

**Protocol for dsDNA recombineering to insert a drug cassette into the chromosome
or a BAC**

- 1.** Electroporate 100-300ng of purified drug cassette-containing PCR product into freshly prepared electrocompetent, recombineering-proficient cells.
 - 2.** After the 1ml L broth is added and the mix is transferred to a sterile culture tube, outgrow the cell for 2-3 hours at 32° with rolling or shaking. This allows for expression of the drug cassette as well as segregation of chromosomes (or BACs) without the cassette.
 - 3.** After outgrowth, make 10-fold serial dilutions in TMG (recipe below), minimal salts or similar osmotically balanced medium. Plate 0.1ml of dilutions on L plates containing the appropriate drug (see FAQ for drug concentrations to use). As you expect approximately 10^3 - 10^4 recombinants per 10^8 surviving cells, then either the undiluted, 10^{-1} or 10^{-2} dilutions should yield a good number of drug resistant recombinants. Plates should be incubated at 32° to prevent further expression of the Red functions.
- Note:** Total cells that survive electroporation can be determined by plating dilutions on L plates. Expect approximately 10^8 /ml although we have seen up to 10-fold fewer with some recombineering strains.
- 4.** Several recombinants should be struck out for single colonies on the appropriate L+drug plates.
 - 5.** The drug cassette insertion should be confirmed by PCR across both new junctions (see FAQ). Make certain there is not also bands for the WT gene, indicating a duplication.

TMG

10 mM Tris base
10 mM MgSO₄
0.01% gelatin