

## Preparing Cells Competent for Recombineering

*Overview* The preparation of recombineering-proficient cells that are ready for electrotransformation

*Duration* About 3 hours

*Preparation* The previous day, grow a 5 ml overnight culture of the chosen recombineering cells at 30-32°C. Include the appropriate drug if a plasmid is supplying the Red functions.

**Caution** *Do not grow recombineering cells at temperatures greater than 34°C.*

**Caution** *Maintain sterile technique throughout the rest of the protocol.*

**1** Dilute the overnight culture by adding 0.5 ml of the overnight to 35 ml of LB medium with the appropriate drug(s) if needed, in a 250 ml (or 125 ml) baffled Erlenmeyer flask. Dilute the overnight at least 70-fold. Grow cells in a H<sub>2</sub>O bath at 32°C with shaking (200rpm) until the OD<sub>600</sub> is from 0.4-0.5 (approximately 2 hrs).

*Tip* *Cells with different genotypes will grow at different rates. Having the proper OD<sub>600</sub> is critical – the recombination will not work if the density is too high.*

*Tip* *Only add drug to the LB if it is needed to maintain a plasmid.*

**2** Transfer half the culture to a 50 ml baffled Erlenmeyer flask and place that flask in a 42°C H<sub>2</sub>O bath to shake at 200rpm; keep the other flask at 32°C. Shake for 15 min. The culture at 42°C is now induced for the recombination functions and the 32°C culture is the uninduced control. Both flasks will be processed identically during the rest of the protocol.

*Tip* *If you do not have baffled flasks, use a 125 ml or larger flask.*

**3** Immediately after the 15 min induction, rapidly chill both cultures in an ice-water slurry; swirl the flasks gently. Leave on ice for 5-10 min. Label and chill the necessary number of 35-50 ml centrifuge tubes for the induced and uninduced cells.

*Tip* *Pre-chill the sterile distilled H<sub>2</sub>O that will be used for washes. Keep 200 ml bottles of distilled water at 4°C for this purpose and put it on ice as needed. Also chill electrotransformation cuvettes and microcentrifuge tubes for later parts of this step.*

**4** Transfer both the induced and uninduced cultures to the chilled centrifuge tubes and centrifuge 7 min at ~6500 x g (6700 rpm in a Sorvall SA-600 rotor) at 4°C. Using sterile technique, aspirate or pour off supernatant.

**5** Add 1 ml ice-cold sterile distilled H<sub>2</sub>O to the cell pellet and gently suspend cells with a large disposable pipet tip (do not vortex). After cells are well suspended, add another 30 ml of ice-cold distilled H<sub>2</sub>O to each tube, seal, and gently invert to mix, again without vortexing. Centrifuge tubes again as in previous step.

**6** Promptly decant the 30 ml supernatant very carefully from the soft pellet in each tube and gently suspend each cell pellet in 1 ml ice-cold distilled H<sub>2</sub>O.

*Tip* As the pellets are very soft, tubes must be removed promptly after centrifugation and care should be taken not to dislodge the pellet. It is ok at this step to leave a small amount of supernatant in the tube.

**7** Transfer the suspended cells to pre-chilled microcentrifuge tubes. Centrifuge 30 sec at maximum speed in a 4°C refrigerated microcentrifuge. Carefully remove the supernatant and suspend cells in 200 µl sterile ice-cold distilled H<sub>2</sub>O and keep on ice until used.

*Tip* This protocol will prepare enough cells for four or five electroporations. If more cells are needed, the best way to do so is to prepare additional flasks.